

NEW Life/form® Multi-Venous IV & Injection Arm

Product Name	NEW Life/form® Multi-Venous IV & Injection Arm
Addresses Training Needs/Curricula	Addresses training requirements for certification or graduation: <ul style="list-style-type: none"> • Lab Technician • Licensed Practical Nurse • Registered Nurse • Nurse Practitioner • Anesthesiology • EMT, Paramedic and advanced degrees Suitable for long-term care facilities and emergency room training.
Configuration	<ul style="list-style-type: none"> • 2 skin tones (light or medium). • Durable and long-lasting with replaceable skin and veins. • Noticeable "pop" as needle enters vein. • Self-sealing veins and skin. • Complex vein structure. • Easy to clean and maintain.
Product Code/SKU	LF01271 - Light skin tone LF01271A - Light skin tone w/ continuous circulation pump LF01272 - Medium skin tone LF01272A - Medium skin tone w/ continuous circulation pump
Accessories included	3 cc syringe with needle* 12 cc syringe with needle* 2 fluid supply bags Infusion butterfly needle* Pinch clamp 2 small towels Infusion needle Blood powder bottle (2 pints) *Needles included are for training purposes only and do not meet standards for use on actual patients.
Replaceable items	Replacement skin (light) - LF01274 Replacement skin (medium) - LF01275 1 qt. simulated blood - LF00845 1 gal. simulated blood - LF00846 Shoulder injection pad replacement - LF00991 Fluid supply stand - LF01022 Fluid supply bag - LF01130 Replacement tubing kit - LF01277 Nasco cleaner - LF09919
Dimensions	Arm: 30" long x 6-1/2" wide x 4-1/2" deep Weight without packaging: 5.90 lbs Packed ship weight: 17 lbs

Be READY



Build critical skills for when it counts



5-year Warranty

NEW Life/form® Multi-Venous IV & Injection Arm

A new and improved evolution of our sophisticated Life/form® Venipuncture and Injection Arm with complete venous access. This trainer simulates a complete full-size adult male arm.


- Upgraded with soft, lifelike skin “feel” that is long-lasting and durable.
- Narrowed floating vein site.
- Intramuscular and intradermal injections.
- Blood drawing and intravenous infusion.
- Primary vein, secondary vein and IV access.
- Non-visible/non-palpable vein access.
- Rolling, disappearing and “difficult” vein access simulation.




Replaceable skin and veins for long-term usability



Noticeable “pop” as needle enters vein



8 fluid lines



Deltoid muscle injection sites





Soft, flexible fingers and wrist

Available in light or medium skin tones



Realistic skin texture



Complete venous access: Hand veins

- Dorsal metacarpal veins
- Dorsal Intercapitular veins
- Branch of cephalic vein of thumb

Arm veins

- Basilic vein
- Cephalic vein
- Median cephalic vein
- Median basilic vein
- Median cubital vein

- Median antebrachial vein
- Accessory cephalic vein

Optional continuous circulation pump

This product is intended to be used as an instructional aide for teaching and demonstrating IV and injection procedures. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies or function as needed.