

MW48

Full Perinatal Patient Care Simulator "Hana"

Be a part of Hana's journey as a new mother

Training manikin from gestational to post-natal period
 Let us help you take your first step in perinatal care!



FEATURES

- 1 | Interchangeable abdomen to simulate perinatal examination, ultrasound examination of 23-week-fetus, cervical examination, labor, perineal tear and puerperal.

POINTS




One-touch interchangeable parts



Compatible with birthing chairs




Durable yet realistically elastic



Accurate anatomy

SKILLS



Perinatal examination

During pregnancy

Perinatal examination

- | Palpation (e.g. using Leopold's maneuvers, cephalopelvic disproportion assessment)
- | Repositionable fetus
- | Auscultation of fetal sounds, with adjustable heart rates from 60 to 180
- | Measurable abdomen girth, fundal height and external pelvimetry



Ultrasound examination

- | Screening of 23-week-fetus
- | Fetal size assessment using BPD, AC, FL, and so on
- | Measurement of amniotic fluid volume
- | Determination of fetus presentation



Childbirth

Cervical examination

- | 5 types of cervical dilation models
- | Determination of bishop score (cervix score)
- | Check fetal head position



During labor

- | Various positions: Supine, lateral, all-fours
- | Perineal protection
- | Cephalic and breech presentations
- | Obstetric forceps delivery and vacuum extraction delivery
- | Tying and cutting of the umbilical cord
- | Delivery of the placenta
- | Inspection of the amniotic sac
- | Insertion of urinary catheter



Perineal tear

- | Perineorrhaphy:
 - single interrupted suture
 - vertical mattress suture
- | Repetitive suture practice for first degree perineal tear



Puerperal

Puerperal assessment

- | Assessment uterine involution in early puerperium: 4 types of variations
- | Measurement and palpation of uterine fundus
- | Perineal cleaning and assessment
- | Postpartum massage
- | Assessment of anal prolapse



DESCRIPTIONS

SPECIFICATIONS

Size: 164 cm/ 63 inch Weight: 20 kg / 2.2 lbs

SET INCLUDES

* Components for basic manikin. Training modules are to be added according to requirements

1 manikin	1 injection pad (median vein)
1 wig	1 irrigation bag
1 face mask for Hana	1 drainage pump /drainage hose
1 denture	1 lubricant for training models
1 female genitalia	1 talcum powder
1 maternity clothes	instruction manual

MATERIALS

Soft resin
Latex free

REPLACEMENT PARTS

11229-050	lubricant for training models	11406-040	wig
11251-030	2 shoulder IM pads	11415-030	genitalia unit
11416-040	5 velamen sheets	11415-010	cervical dilation inserts
11416-030	10 umbilical cords for omphalotomy	11416-020	placenta
11417-010	2 genitalia insert		
11415-040	lubricant (500ml)		

RECOMMENDED DEVICES

-