



PATIENT CARE



CAREGIVER SUPPORT



BED ACCESSORIES AT A GLANCE



OUR BED ACCESSORIES RANGE IS DESIGNED TO MEET THE NEEDS OF CARE HOME AND HOSPICE RESIDENTS AND THOSE BEING CARED FOR AT HOME. ALL SOLUTIONS CAN WORK IN CONJUNCTION WITH YOUR SPECIFIC NEEDS, REQUIREMENTS AND RISK ASSESSMENTS.

THIS AT A GLANCE GUIDE AIMS TO PROVIDE YOU ALL THE KEY INFORMATION YOU NEED TO SUPPORT YOUR PRODUCT CHOICE, ALL IN ONE PLACE.



Bed Rail Protector Range

Designed to be fixed in place to support residents who may be prone to tossing and turning when in bed, often injuring themselves in the process. Most profiling beds have two full length rails, and our Bed Rail Protectors Range has been designed to fit over these rails.



MRSA Resistant - with the use of a non-porous PVC outer fabric



FR Tested - Flame Retardant tested to British Standard BS 7177



Wipe Clean Finish - supports infection control procedures and is easy to maintain



White Colour - to promote a fresh looking, contemporary environment

Bed Rail Protectors



These durable bed rail protectors are supplied as a pair and provide a cushioned barrier for residents

- Available in standard and full length versions
- White PVC with foam filling
- Hook and loop fastening, for ease of fitting and maintenance

Mesh Bed Rail Protector



As well as providing a padded barrier, the mesh inserts improve visibility between each side rail to reduce the feeling of entrapment

- Available in full length
- White PVC with foam filling with mesh inserts
- Zipper fastening, for durability

Connected Bed Rail Protector



The mesh base helps to secure the protector sides in place, forming a safety net to help prevent entrapment between the mattress and the bottom of the side rails / protector sides

- Available in full length
- White PVC with foam filling and mesh connecting sheet
- Multi-fitting mesh base allows for the movement of a profiling bed
- Hook and loop fastening, for ease of fitting and maintenance



Bed Support Accessories Range

Made from blue PVC with a high density foam filling, our range of bed support accessories are designed to be robust and durable to ensure peace of mind for users and carers alike. All Products in this range are FR tested and can be wiped clean.



FR Tested - Flame Retardant tested to British Standard BS 7177



Wipe Clean Finish - supports infection control procedures and is easy to maintain

Side Wedges With Connecting Sheet



The high density foam filled pair of side wedges with a breathable connecting sheet are designed to minimise the risk of residents falling out of beds where rails are not used

- Blue PVC with high density foam filling
- Hook and loop fastening of the wedges to connecting sheet
- Durable, quick release straps for a secure fit to the mattress
- Wedges and connecting sheet also available separately

Multi-Purpose Bed Wedge



A comfortable, versatile positioning aid that can be used upright or horizontally to encourage a healthy postural position when sitting or reclining in bed. Alternatively it can be used for leg support. Further positioning options can be created using two bed wedges

- Flexible usage
- Blue PVC with high density foam filling


Fall-Out Mat



Designed for use with low-level beds to reduce the risk of injury from falling out of bed. Our Fall-Out Mat features durable carry straps and folds away for ease of storage

- Blue PVC with high density foam filling
- Durable carry straps
- Folds away for easy storage

Static Mattress Extension Infill



This high density foam infill is designed to fill the gap often left at the foot of extended profiling beds, enhancing comfort for taller residents

- Blue PVC with high density foam filling

MIP BED ACCESSORIES

PRODUCT SPECIFICATIONS

Bed Rail Protectors

Product Code	Dimensions	Size
135/W/PU	137 x 86 x 2cm	Standard Length
136/W/PU	198 x 86 x 2cm	Full Length

Connected Bed Rail Protector

Product Code	Dimensions	Size
143/W/PU	198 x 86 x 2cm	Full Length

Mesh Bed Rail Protector

Product Code	Dimensions	Size
144/W/PU	201 x 86 x 2cm	Full Length

Side Wedges with Connecting Sheet

Product Code	Dimensions	Product
138/W/PU	90 x 20cm	Side Wedges with Connecting Sheet
138/WEDGES	90 x 20cm	Side Wedges
137/W/PU	90 x 90cm	Connecting Sheet

Multi-Purpose Bed Wedge

Product Code	Dimensions
141/W/PU	60 x 24 x 57cm

Fall Out Mat

Product Code	Dimensions
139/W/PU	190 x 60 x 5cm

Static Mattress Extension Infill

Product Code	Dimensions
140/W/PU	90 x 20 x 15cm