



# Product Specification

## Multi-adapter, Luer



### Product description

Order number	14.1205
Product description	Multi-adapter, Luer, connection: Luer, for S-Monovette® , 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free
PZN number	6796193

### Product characteristics

Design	Luer
Suitable for	S-Monovette®
Connection	Luer

### Size

Length of product	33 mm
-------------------	-------

### Material & colours

Product material	Polyethylene therephthalate (PET)
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-02-07

This document was prepared by EDP support and is valid without signature.

## Product Specification Multi-adapter, Luer

### Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of box	268 mm
Width of box	81 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	178 mm
Height of case	287 mm
Case volume	0.0191 cbm
Weight of product	0.0014 kg
Weight of case	2.34 kg
EAN of inner box	4038917085793
EAN case	4038917004428

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.