


MARSDEN

MEDGUARD[®]
Professional Healthcare Supplies

V-150 Large Veterinary Scale



- Capacity: 150kg or 300kg
- Accurate from: 50g or 100g
- Stainless steel weighing scale
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty

Marsden V-150:

Large Veterinary Scale

- Capacity to 150kg or 300kg
- Accurate from 50g
- Portable scale for weighing pets
- Lightweight and easy to transport
- Automatic hold for accuracy
- Perfect size for cats and dogs
- Weigh a wide range of animals with ease
- Wipe down top
- Easy to clean and keep hygienic
- Simple and easy to use
- Tare to remove unwanted weight
- Covered by 1 year warranty



What is the V-150?

The V-150 is a large, professional veterinary scale which has a large stainless steel base. It is accurate to 50g or 100g and comes as standard with a non-slip mat.

The V-150 has quickly become a best-selling animal weighing scales to veterinary surgeries as it is accurate enough to weigh many different animals and you will often find them situated in the main reception.

Easy to clean

To ensure that longevity of this animal weighing scale, we designed this model to be have a stainless steel construction, meaning it can be cleaned easily.

This means that animals can stand on the scale when they are wet and muddy, with no detrimental effect to the weighing scale and it can be wiped down by the veterinary staff and it is good to go.

Highly portable

We understand that whether the scale is been used at home or within a veterinary surgery, it is highly likely that the scale will need to be transported and moved from time to time.

For this reason, the V-150 comes as standard with in-build handles and wheels to ensure that the user can easily manoeuvre the scale into position.

Capacity:	150kg or 300kg	Power Supply:	Mains/Rechargeable
Graduations:	50g or 100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	900mm x 550mm x 80mm	Optional Extras:	N/A
Weight of Scale:	15kg	Warranty Period:	1 Year