

Sani-Cloth® 70% is an ideal disinfectant for use in fast-paced environments requiring a short contact time and broad efficacy coverage against microorganisms.

Tuberculocidal, virucidal, bactericidal, and yeasticidal in 2 minutes.



Recommended Uses:

- + General surfaces
- + Common equipment & non-invasive medical devices
- + As part of a two-step process with Sani-Cloth® Detergent to clean prior to disinfection.

+ Ready to Use

Pre-dosed, low linting disposable alcohol disinfection wipes, requiring no mixing or measuring, with no smear formulation.

+ Maximum Flexibility

Can be used on non-invasive medical devices; non-porous hard surfaces.

+ Supporting Resources

Point of Use accessories including canister bracket holders, "I Am Clean" tape & sticky notes, IFU videos and posters.

+ Proven Efficacy

Complete disinfection in 2 minutes including Coronavirus and Flu in 15 seconds.

+ Fast Acting

Quick drying, no residue formulation containing 70% Isopropanol alcohol, allowing the disinfected area or equipment to be used shortly after application.

+ Supports Carbon Net Zero Goals

Manufactured in the UK, resulting in supply chain reliability and reducing carbon footprint.

Product Efficacy

	Test Method	Organism	Contact Time	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	2 min (d)
	EN 16615	Staphylococcus aureus	60 sec (c)	2 min (d)
	EN 16615	Enterococcus hirae	60 sec (c)	
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeast	EN 16615	Candida albicans	60 sec (c)	2 min (d)
	EN 13624	Candida albicans	30 sec (c)	
Enveloped virus	EN 14476	Vaccinia Virus	30 sec (c)	60 sec (d)
	EN 14476	Corona,HIV, Hep & C, SARS, Flu	15 sec (c)	
TB	EN 14348	Mycobacterium terrae	60 sec (c)	

*Clean Conditions (c) Dirty Conditions (d)

Product Details



Format	Wipe Size	Wipe Count	QTY per Case	Classification	Product Code
200ct Large Canister	200mm x 200mm	200	6	MD Class IIa	E701020