



Ascender Telescopic ROM knee Brace

Product Code: ASC/10

Designed and manufactured in the UK by Promedics. The new Ascender telescopic ROM knee brace is universal, easily adjustable and secure. With additional features offering superior comfort and support to patients recovering from knee surgery, injuries or instabilities.

- Universal sizing
- Maximum strength, minimum weight
- High grade anodised aluminium
- Length range 44 cm – 70 cm
- 5 drop lock settings
- Quick-clip Range of Motion Settings

The range of motion settings can be easily adjusted via the spring loaded quick-clips by a simple 'pull and slide' motion.

The design prevents unintentional changes to the flexion and extension settings.

The drop lock is ergonomically designed to be easily locked and unlocked due to its curved concave profile.

Size	Product Code
Universal	ASC/10

Ascender C/W Foam Pads

Product Code: ASCFP/10

Size	Product Code
Universal	ASCFP/10



Scan QR Code
to watch video

www.ascender.pro

 **ascender**[®]

Range of Motion Control



Universal Sizing
Length range 44cm – 70cm



Designed for Comfort
Optimised for Airflow
Silicone for Ultra Grip

