

Jura Black Wrist Thumb Brace

The Jura Black Wrist Thumb Brace immobilises and protects the wrist following tendonitis and other soft tissue injuries such as sprains/ strains, and post-trauma healing.

Lightweight, breathable and contoured for a comfortable fit, with an adjustable palmar bar. Slide-on design with D-ring closures allow easy application.



Durable, tough webbing bar casing and a two-piece brace with longitudinal seam conforms perfectly to the wrist.



Slide-on design with D-ring closure allows for easy application.



Made from dual compression elastic for a great fit, the Neoprene fabric is breathable and comfortable against the skin.

Size	Left	Right	Wrist Circ.
XS	JRWTB/12	JRWTB/22	< 15cm
S	JRWTB/13	JRWTB/23	15cm - 17cm
M	JRWTB/15	JRWTB/25	17cm - 19cm
L	JRWTB/17	JRWTB/28	19cm - 22cm
XL	JRWTB/18	JRWTB/28	23cm - 25cm