

Jura Black Wrist Brace - Regular / Long

The Jura Black Wrist Brace immobilises and protects the wrist following tendonitis and other soft tissue injuries such as sprains/ strains, and post-trauma healing.

Lightweight, breathable and contoured for a comfortable fit, with an adjustable palmar bar. Slide-on design with D-ring closures allow easy application.



Durable, tough webbing bar casing and a two-piece brace with longitudinal seam conforms perfectly to the wrist.



Slide-on design with D-ring closure allows for easy application.



Made from dual compression elastic for a great fit, the Neoprene fabric is breathable and comfortable against the skin.

Size	Regular - Left	Regular - Right	Long - Left	Long - Right	Wrist Circ.
XS	JRWB/12	JRWB/22	JRWBL/12	JRWBL/22	< 15cm
S	JRWB/13	JRWB/23	JRWBL/13	JRWBL/23	15cm - 17cm
M	JRWB/15	JRWB/25	JRWBL/15	JRWBL/25	17cm - 19cm
L	JRWB/17	JRWB/28	JRWBL/17	JRWBL/27	19cm - 22cm
XL	JRWB/18	JRWB/28	JRWBL/18	JRWBL/28	23cm - 25cm