

SK Knee

The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity. Hyperstop hinges provide medial and Lateral support for the unstable knee, while allowing flexion.

- + Mild instabilities of the knee
- + Meniscus injuries
- + Trauma or degenerative deformities of the knee joint
- + Protection in the later stages of rehab

Available in black or beige versions.



Size	Beige	Black	Knee Circ.
XS	SK/02	SK/02/BK	25cm - 28cm
S	SK/03	SK/03/BK	28cm - 34cm
M	SK/05	SK/05/BK	34cm - 40cm
L	SK/07	SK/07/BK	40cm - 46cm
XL	SK/08	SK/08/BK	46cm - 56cm
XXL	SK/09	SK/09/BK	56cm - 61cm