



Upper arm



Model	EVOLV	M7 Intelli IT	M6 Comfort	M4 Intelli IT	M3 Comfort	M3	M2	M2 Basic
Product Code	HEM-7600T-E	HEM-7361T-EBK	HEM-7360-E	HEM-7155T-EBK	HEM-7155-E	HEM-7154-E	HEM-7121-E	HEM-7120-E
All in One No tubes. No wires. Simply accurate.	●							
Intelli Wrap Cuff Technology Accurate results in any position around the upper arm	●	●	●	●	●			
OMRON connect	●	●		●				
Intellisense Technology Measurement without unnecessary over-inflation of the cuff	●	●	●	●	●	●	●	●
Atrial Fibrillation Indicator		●	●					
3x Measurement Function		●	●					
Comparison Function		●	●					
Irregular Heartbeat Detection	●	●	●	●	●	●	●	
Blood Pressure Indicator		●	●	●	●	●		
Blood Pressure Level Indicator Indicates if result is above recommended level							●	
Body Movement Detection Signals if excessive movement detected during measurement	●	●	●	●	●	●		
Cuff Wrap Guide Supports accurate measurement	●	●	●	●	●	●	●	
Cuff - Type & Size	Intelli Wrap Cuff (Preformed) 22-42 cm	Intelli Wrap Cuff (Preformed) 22-42 cm	Intelli Wrap Cuff (Preformed) 22-42 cm	Intelli Wrap Cuff (Preformed) 22-42 cm	Intelli Wrap Cuff (Preformed) 22-42 cm	Easy Cuff 22-42cm	Medium Cuff 22-32 cm	Medium Cuff 22-32 cm
Optional Cuff						Small Arm Cuff 17-22cm	Small Arm Cuff and Easy Cuff 17-22cm and 22-42cm	Small Arm Cuff and Easy Cuff 17-22cm and 22-42cm
Memory Capacity	1 user x 100 *	2 users x 100	2 users x 100	2 users x 60	2 users x 60	2 users x 60	1 user x 30	1 user x 1
Storage Case Included	yes	yes	yes	yes	yes	yes		

* This device stores 100 measurements in its memory. To view your measurements you need to synchronize the device with your smart phone using the OMRON connect app.