




Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Paediatric and Infant Stethoscopes	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.		Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.
Approximate Length	71cm	71cm	69cm	69cm	69cm	69cm
Approximate Weight	118g	105g, 95g	150g	177g	185g	232g
Warranty	2 years	3 years	5 years	7 years	7 years	2 years
Diaphragm	Tunable	Non-Tunable (Floating)	Tunable (2)	Tunable (2)	Tunable	Tunable (2)
Approximate Chestpiece Diameter	Double-sided 5.2 x 4.5cm 3.5 x 3.2cm	Double-sided Paediatric 3.4cm 2.7cm Infant 2.8cm 1.7cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Single-sided 4.8cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm
 Features	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital Connects to Eko software on a smart device** to capture, store and share data

The power to learn,
in the palm of your hand.

3M™ Littmann® Learning App – Download free version in the App Store® or Google Play™ Store
 Need more help? Visit Littmann.co.uk

*At peak frequency, vs. analogue mode
 **Not included

3M United Kingdom PLC
 Charnwood Campus | 10 Bakewell Road | Loughborough | Leicestershire LE11 5RB