

Technical Data

Evacuation Sledge

TM Ref:

EVSG-EC

Description

This device consists of an inner foam layer for patient comfort. There is strong webbing and a foot pocket to hold the patient in place, pulling loops at each end and a tough slippery coated bottom side. The loop at the head end is to hold back when descending stairways. The whole item is stored in a pouch and can hang on the wall. One or two people required.

Dimensions

Size 1900 x 600 x 40 mm

Foam Inner

Foam comfort layer

Materials Used

Webbing 50 mm in width

Buckles

50 mm Side Release Buckle (plastic)

Fabrics

Woven Poly / Cotton

Base Fabric

High tenacity Polyester

Care instructions

For general cleaning purposes, wipe the cover with warm water containing a neutral detergent. Where disinfection is required following contamination, refer to hospital disinfection policy. Do not use phenolic based cleaners or bleach on this product .If laundering is necessary and disinfection is required wash up to at 40°C. Tumble dry at low temperature. Not iron.



Wash 40°, Do not bleach, Tumble dry on low temp, Do not iron