

Name	Value
Undercounter pharmacy refrigerator	MD14DC-P301QX23XXXZ1SN4B
EAN code	4770104014055
Gross capacity, l	130
Net capacity, l	100
Factory temperature setting, °C	+5
Fluctuation, °C	0,5
Gradient, °C	2,9
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.
Recovery time after 1 min of door opening (empty, +25 °C TU, according to EN 60068-3)	8 min.
Positive temperature deviation from the setpoint	0,5
Negative temperature deviation from the setpoint	-1,4
Heat output, kJ/h	135
Ambient temperature range, °C	from +10 to +35
Airborne acoustical noise emissions, dB(A) re1 pW	42
Refrigerant	R600a
Daily energy consumption, kWh / 24 h	0,751
Annual energy consumption	274
Rated current, A	0,5
Rated voltage, V	220 - 240
Rated frequency, Hz	50
Cooling system	Fan forced cooling
Defrost	Automatic
Volt-free alarm contact	No
Self closing door	Yes
Access port	Yes
Interior light	LED ceiling light strip
Malfunction warning signal	Optical and acoustical
Protection against low temperature	Yes
Number of shelves	2 adjustable, 2 bottom
Shelf size (W x D), mm	420 x 324, bottom: 465 x 263 and 445 x 180
Maximum load per shelf, kg	20
Side wall material	Steel, white painted
Inner cabinet material	HIPS polystyrene
Type of door	Double panel selective glass in aluminium profile

Name	Value
Type of lock	Mechanical
Adjustable foot at the front	2
Foot behind	2
External dimensions (H x W x D), mm	870 x 560 x 620
Internal dimensions (H x W x D), mm	725 x 465 x 392
Package dimensions (H x W x D), mm	950 x 570 x 640
Weight net, kg	37
Weight gross, kg	40
Pharmacy refrigerator packed on a pallet	MD14DC-P301QX23XXXXZ1SN4A
EAN code	4770104014055
Dimensions (H x W x D), mm	1080 x 570 x710
Weight gross, kg	48,5