

3M™ Littmann® Stethoscopes

Introducing the next
generation of a clinical icon.



Open bell versatility.

The pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

Introducing the

3M™ Littmann® Classic III™ Stethoscope

Right out of the box, it lives up to its name.

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The new Classic III stethoscope offers you:



High acoustic sensitivity for performing general physical assessments.



Tunable diaphragms on both the adult and pediatric sides of the chestpiece.



Single-piece diaphragm that's easy to attach and easy to clean because the surface is smooth without crevices.



Next-generation tubing with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.

Here's what makes the
3M™ Littmann® Classic III™ Stethoscope stand out.

Features

Tunable technology is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

Soft-sealing eartips give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

Open bell stays clear of dirt and debris by keeping it covered with the small diaphragm.

Adjustable headset lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

No natural rubber latex or phthalate plasticizers are used in the tubing or any other component—to help protect human health and the environment.

Made in the USA for consistent quality. We apply rigorous process and quality controls to every Littmann stethoscope. You'll always get a durable, dependable, clinical tool that will serve you for years to come.

Classic III stethoscope specifications

Length	69cm / 27.2 inches
Weight	150g
Acoustic performance	7
Warranty	5 years
CHESTPIECE	
Construction	Two-sided, adult and pediatric
Material	Stainless steel
Diaphragm type	Tunable, single-piece
Diaphragm diameter	Adult: 4.3cm / 1.7 in Pediatric: 3.3cm / 1.3 in
TUBING	
Construction	Single lumen (single acoustic channel within tube)
Material	Not made with natural rubber latex or phthalate plasticizers
HEADSET	
Material	Wide diameter aerospace alloy / anodized aluminum

Colors and Finishes

We periodically introduce new tubing colors and chestpiece finish combinations for the Classic III stethoscope. Visit littmann.com for the latest product information.

Due to limitations of the printing process, the colors shown are approximations of the actual tubing colors.



Available tubing colors and finish combinations

TUBING COLOR	STANDARD FINISH	SPECIAL FINISH
Black	5620	Black 5803
Navy Blue	5622	
Raspberry	5626	
Burgundy	5627	
Chocolate		Copper 5809
Lavender	5832	
Turquoise	5835	
Lemon-Lime	5839	
Emerald	5840	
Rose Pink*	5639	

*Not available in the U.S.